

Anti-HAPLN1 Antibody Picoband® PE Conjugated

Catalog Number: A05980-3-PE

About HAPLN1

Predicted to enable hyaluronic acid binding activity. Predicted to be an extracellular matrix structural constituent conferring compression resistance. Predicted to be involved in nervous system development; positive regulation of neuroblast proliferation; and skeletal system development. Located in collagen-containing extracellular matrix.

Overview

Product Name	Anti-HAPLN1 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P10915

Technical Details

Immunogen	E.coli-derived human HAPLN1 recombinant protein (Position: D16-N354). Human HAPLN1 shares 96.5% amino acid (aa) sequence identity with both mouse and rat HAPLN1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HAPLN1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.