

## Anti-VAMP1 Antibody

Catalog Number: A05982-2


### Overview

Product Name	Anti-VAMP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-VAMP1 Antibody catalog # A05982-2. Tested in WB,IHC,ICC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23763


### Technical Details

Immunogen	Full length recombinant protein of human VAMP1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC/IF: 1:50-1:200 IHC: 1:50-1:200

## Anti-VAMP1 Antibody (A05982-2) Images



Western blot analysis of VAMP1 on different tissue lysates using anti-VAMP1 antibody at 1/2,000 dilution. Positive control: Lane 1: Mouse cerebellum Lane 2: Rat brain



ICC staining VAMP1 in 293T cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-VAMP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.