

Anti-VAMP1 Antibody

Catalog Number: A05982-3

About VAMP1

VAMP1, also named as synaptobrevin 1, is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. VAMP1 is a part of the SNARE (soluble NSF-attachment protein receptor) complex. Characterized by a common sequence called the SNARE motif, SNARE proteins are involved in membrane fusion and vesicular transport (PMID: 11252968). VAMP1 is involved in synaptic vesicle exocytosis, a fundamental step in neurotransmitter release.

Overview

Product Name	Anti-VAMP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-VAMP1 Antibody catalog # A05982-3. Tested in WB, IHC, ICC/IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P23763

Technical Details

Immunogen	E.coli-derived human VAMP1 recombinant protein (Position: A3-M97).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-VAMP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.