

Anti-TMEPA PMEPA1 Antibody

Catalog Number: A05994

About PMEPA1

Involved in down-regulation of the androgen receptor (AR). Enhances ubiquitination and proteasome-mediated degradation of AR, probably by recruiting NEDD4.

Xu L.L., Genomics 66:257-263(2000).

Rae F.K., Mol. Carcinog. 32:44-53(2001).

Deloukas P., Nature 414:865-871(2001).

Overview


Product Name	Anti-TMEPA PMEPA1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-TMEPA PMEPA1 Antibody (Catalog # A05994). Tested in WB, IHC, IF applications. This antibody reacts with Human.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q969W9

Technical Details


Immunogen	Synthesized peptide derived from internal of human TMEPA.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	Western blotting: 1:500~1:3000

Immunohistochemistry: 1:50~1:100
Immunofluorescence: 1:100~1:500


Anti-TMEPA PMEPA1 Antibody (A05994) Images



Immunohistochemistry analysis of paraffin-embedded human prostate tissue using TMEPA antibody A05994.



Western blot analysis of extracts from HT-29 cells, using TMEPA antibody A05994.



Immunofluorescence analysis of HUVEC cells, using TMEPA antibody A05994.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TMEPA PMEPA1 Antibody