

## Anti-Adenylate cyclase type 9 ADCY9 Antibody

Catalog Number: A06038-1


### Overview

Product Name	Anti-Adenylate cyclase type 9 ADCY9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Adenylate cyclase type 9 ADCY9 Antibody catalog # A06038-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60503


### Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human ADCY9.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200


## Anti-Adenylate cyclase type 9 ADCY9 Antibody (A06038-1) Images



Western blot (WB) analysis of ADCY9 polyclonal antibody at 1:500 dilution Lane1:Hela cell lysate Lane2:sp2/0 cell lysate Lane3:PC12 cell lysate



Immunohistochemistry (IHC) analyzes of ADCY9 pAb in paraffin-embedded human breast carcinoma tissue at 1:50, showing cytoplasm and membrane staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.



Immunohistochemistry (IHC) analyzes of ADCY9 pAb in paraffin-embedded human stomach carcinoma tissue at 1:50, showing cytoplasm and membrane staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

BS5590  
Lot: CN75131

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Adenylate cyclase type 9 ADCY9 Antibody

For Research Use Only. Not for use in diagnostic procedures.