

Anti-eIF3 alpha EIF3J Antibody

Catalog Number: A06040

Overview


Product Name	Anti-eIF3 alpha EIF3J Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-eIF3 alpha EIF3J Antibody catalog # A06040. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75822

Technical Details

Immunogen	Synthesized peptide derived from eIF3alpha . at AA range: 40-120
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:40000


Anti-eIF3 alpha EIF3J Antibody (A06040) Images

SH-SY5Y 293T 3T3



Tubulin 55KD

Western Blot analysis of various cells using Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000



eIF3α 30KD

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-eIF3 alpha EIF3J Antibody

For Research Use Only. Not for use in diagnostic procedures.