

Anti-SR-1D HTR1D Antibody

Catalog Number: A06046


Overview

Product Name	Anti-SR-1D HTR1D Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SR-1D HTR1D Antibody catalog # A06046. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P28221


Technical Details

Immunogen	Synthesized peptide derived from SR-1D . at AA range: 120-200
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using a epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 ELISA 1:10000 IF 1:50-200

Anti-SR-1D HTR1D Antibody (A06046) Images



Western Blot (WB) analysis of HepG2 823 293T AD293 HeLa using SR-1D antibody.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Brain, antibody was diluted at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SR-1D HTR1D Antibody

For Research Use Only. Not for use in diagnostic procedures.