

Anti-Kallikrein-11 KLK11 Antibody

Catalog Number: A06052-1


Overview

Product Name	Anti-Kallikrein-11 KLK11 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Kallikrein-11 KLK11 Antibody catalog # A06052-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UBX7


Technical Details

Immunogen	Recombinant fusion protein of human KLK11(NP_001129504.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:100-1:200

Anti-Kallikrein-11 KLK11 Antibody (A06052-1) Images



Immunohistochemistry of paraffin-embedded mouse brain using KLK11 antibody at dilution of 1:100 .Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded mouse lung using KLK11 antibody at dilution of 1:100 .Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Kallikrein-11 KLK11 Antibody

For Research Use Only. Not for use in diagnostic procedures.