

Anti-SYCP1 Antibody Picoband® Biotin Conjugated

Catalog Number: A06079-1-Biotin

About SYCP1

Enables double-stranded DNA binding activity. Involved in protein homotetramerization. Predicted to be located in synaptonemal complex. Predicted to be active in central element; male germ cell nucleus; and transverse filament.

Overview

Product Name	Anti-SYCP1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q15431

Technical Details

Immunogen	E.coli-derived human SYCP1 recombinant protein (Position: D95-H717). Human SYCP1 shares 75.8% and 76.6% amino acid (aa) sequence identity with mouse and rat SYCP1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SYCP1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.