

Anti-SYCP1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A06079-1-Fluoro488

About SYCP1

Enables double-stranded DNA binding activity. Involved in protein homotetramerization. Predicted to be located in synaptonemal complex. Predicted to be active in central element; male germ cell nucleus; and transverse filament.

Overview

Product Name	Anti-SYCP1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q15431

Technical Details

Immunogen	E.coli-derived human SYCP1 recombinant protein (Position: D95-H717). Human SYCP1 shares 75.8% and 76.6% amino acid (aa) sequence identity with mouse and rat SYCP1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-SYCP1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.