

## Anti-Talin-2 TLN2 Antibody

Catalog Number: A06080

### Overview

Product Name	Anti-Talin-2 TLN2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Talin-2 TLN2 Antibody catalog # A06080. Tested in IHC, IF applications. This antibody reacts with Human, Mouse.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y4G6

### Technical Details

Immunogen	Synthesized peptide derived from part region of human protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	TLN2 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-300 IF 1:50-200

### 1 Publications Citing This Product

1. PubMed ID: 32198023, Baster Z, Li L, Rajfur Z, Huang C. Talin2 mediates secretion and trafficking of matrix metalloproteinase 9 during invadopodium formation. *Biochim Biophys Acta Mol Cell Res.* 2020 Jul;1867(7):118693. doi:10.1016/j.bbamcr.2020.118693. Epub 2020 Mar 18. PMID:32198023

Visit [bosterbio.com/anti-talin-2-tln2-antibody-a06080-boster.html](https://www.bosterbio.com/anti-talin-2-tln2-antibody-a06080-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Talin-2 TLN2 Antibody

For Research Use Only. Not for use in diagnostic procedures.