

## Anti-NUB1 (Y601) Antibody

Catalog Number: A06089Y601


### Overview

Product Name	Anti-NUB1 (Y601) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NUB1 (Y601) Antibody catalog # A06089Y601. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y5A7


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 570-621 of Human NUB1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:10-1:100


## Anti-NUB1 (Y601) Antibody (A06089Y601) Images



Western blot (WB) analysis of NUB1 (Y601) pAb at 1:500 dilution Lane1:HEK293T whole cell lysate(10ug) Lane2:A549 whole cell lysate(10ug) Lane3:U-87MG whole cell lysate(10ug) Lane4:MCF-7 whole cell lysate(10ug) Lane5:PC3 whole cell lysate(10ug) Lane6:PC12 whole cell lysate(40ug) Lane7:AML-12 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of NUB1 (Y601) pAb in paraffin-embedded human tonsil tissue.



Immunoprecipitation of HEK293T cell lysate using NUB1 (Y601) pAb (Sepharose Bead Conjugate) #BD0048(lane 2 and lane 3) .Lane 1 is 30% input. The western blot was probed using NUB1 (Y601) #BS1464. "↑"□supernatant□; "↓(deposition□

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NUB1 (Y601) Antibody

For Research Use Only. Not for use in diagnostic procedures.