

Anti-LMO1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A06093-1-Fluoro594

About LMO1

Rhombotin-1 is a protein that in humans is encoded by the LMO1 gene. This locus encodes a transcriptional regulator that contains two cysteine-rich LIM domains but lacks a DNA-binding domain. LIM domains may play a role in protein interactions; thus the encoded protein may regulate transcription by competitively binding to specific DNA-binding transcription factors. Alterations at this locus have been associated with acute lymphoblastic T-cell leukemia. Chromosomal rearrangements have been observed between this locus and at least two loci, the delta subunit of the T-cell antigen receptor gene and the LIM domain binding 1 gene. Alternatively spliced transcript variants have been described.

Overview

| | |
|----------------------|--|
| Product Name | Anti-LMO1 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P25800 |

Technical Details

| | |
|---------------------|--|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human LMO1, different from the related mouse sequence by two amino acids. |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LMO1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.