

## Anti-SLC27A1 Antibody Picoband® Biotin Conjugated

Catalog Number: A06094-Biotin

### About SLC27A1

Enables biotin transmembrane transporter activity; efflux transmembrane transporter activity; and long-chain fatty acid transmembrane transporter activity. Involved in export across plasma membrane; lipid transport across blood-brain barrier; and monocarboxylic acid transport. Located in plasma membrane.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-SLC27A1 Antibody Picoband® Biotin Conjugated  |
| Reactive Species     | Mouse, Rat   |
| Application          | WB, IHC, ELISA   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.                             |
| Host                 | Rabbit   |
| Uniprot ID           | Q6PCB7   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthetic peptide corresponding to a sequence at the N-terminus of human SLC27A1.   |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | Biotin  |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users.<br>Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.<br>ELISA, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SLC27A1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.