

## Anti-SLC27A1 Antibody Picoband® Cy3 Conjugated

Catalog Number: A06094-Cy3

### About SLC27A1

Enables biotin transmembrane transporter activity; efflux transmembrane transporter activity; and long-chain fatty acid transmembrane transporter activity. Involved in export across plasma membrane; lipid transport across blood-brain barrier; and monocarboxylic acid transport. Located in plasma membrane.

### Overview

Product Name	Anti-SLC27A1 Antibody Picoband® Cy3 Conjugated
Reactive Species	Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6PCB7

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human SLC27A1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-SLC27A1 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.