

Anti-HSPC142 (BABAM1) Rabbit Polyclonal Antibody

Catalog Number: A06132


Overview

Product Name	Anti-HSPC142 (BABAM1) Rabbit Polyclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Rabbit polyclonal anti-BABAM1 Antibody. Catalog# A06132. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Conjugate	Unconjugated
Application	WB
Clonality	Polyclonal
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NWX8

Technical Details

Immunogen	Recombinant protein of human BABAM1
Isotype	IgG
Concentration	2.2mg/ml
Purification	Purified from the immunized serum by affinity chromatography (Protein A/G)
Suggested Dilutions	WB: 1:1000

Anti-HSPC142 (BABAM1) Rabbit Polyclonal Antibody (A06132) Images



Western blot analysis of extracts (35ug) from different cell lines or tissues by using anti-BABAM1 rabbit polyclonal antibody (A06132).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSPC142 (BABAM1) Rabbit Polyclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.