

## Anti-MSRB2 Antibody Picoband® FITC Conjugated

Catalog Number: A06136-1-FITC

### About MSRB2

Predicted to enable actin binding activity; peptide-methionine (R)-S-oxide reductase activity; and zinc ion binding activity. Predicted to be involved in actin filament polymerization and protein repair. Located in mitochondrion.

### Overview

Product Name	Anti-MSRB2 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y3D2

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human MSRB2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-MSRB2 Antibody - FITC**

**For Research Use Only. Not for use in diagnostic procedures.**