

## Anti-DOCK5 Antibody Picoband® HRP Conjugated

Catalog Number: A06138-HRP

### About DOCK5

This gene encodes a member of the dedicator of cytokinesis protein family. Members of this family act as guanine nucleotide exchange factors for small Rho family G proteins. The protein encoded by this gene is thought to associate with adaptors CRK and CRKL, and function in regulation of intestinal epithelial cell spreading and migration on collagen IV. Similar proteins in mouse and zebrafish also function in myoblast fusion.

### Overview

Product Name	Anti-DOCK5 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9H7D0

### Technical Details

Immunogen	E.coli-derived human DOCK5 recombinant protein (Position: Q349-P1776).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DOCK5 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.