

Anti-AMD1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A06146-Fluoro550

About AMD1

S-adenosylmethionine decarboxylase (AdoMet-DC), also known as S-adenosylmethionine decarboxylase proenzyme (SAMDC) or AMD1, is a key enzyme in polyamine biosynthesis. It is localized to chromosome region 6q21-q22. SAMDC has an unusual distribution in polysomes from cells of T lymphocyte origin. It associates predominantly with monosomes and small polysomes with none located in the preribosomal or ribonucleoprotein pool. SAMDC is a critical regulatory enzyme of the polyamine synthetic pathway, and a well-studied drug target. Since SAMDC is a key regulatory enzyme in the synthesis of spermidine and spermine, the marked increase in SAMDC activity in the neonate and the sustained high enzyme levels throughout adulthood, imply a role for these polyamines in both development and mature brain function.

Overview

Product Name	Anti-AMD1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P17707

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human AMD1, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AMD1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.