

Anti-TNFRSF19 Antibody Picoband® PE Conjugated

Catalog Number: A06157-1-PE

About TNFRSF19

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is highly expressed during embryonic development. It has been shown to interact with TRAF family members, and to activate JNK signaling pathway when overexpressed in cells. This receptor is capable of inducing apoptosis by a caspase-independent mechanism, and it is thought to play an essential role in embryonic development. Alternatively spliced transcript variants encoding distinct isoforms have been described.

Overview

| | |
|----------------------|---|
| Product Name | Anti-TNFRSF19 Antibody Picoband® PE Conjugated |
| Reactive Species | Human |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q9NS68 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human TNFRSF19 recombinant protein (Position: E68-L423). Human TNFRSF19 shares 70.6% amino acid (aa) sequence identity with mouse TNFRSF19. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF19 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.