

Anti-POGZ Antibody Picoband® Biotin Conjugated

Catalog Number: A06169-3-Biotin

About POGZ

Pogo transposable element with ZNF domain is a protein that in humans is encoded by the POGZ gene. The protein encoded by this gene appears to be a zinc finger protein containing a transposase domain at the C-terminus. This protein was found to interact with the transcription factor SP1 in a yeast two-hybrid system. Alternatively spliced transcript variants encoding distinct isoforms have been observed.

Overview

Product Name	Anti-POGZ Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q7Z3K3

Technical Details

Immunogen	E.coli-derived human POGZ recombinant protein (Position: D23-R1110). Human POGZ shares 93.7% amino acid (aa) sequence identity with mouse POGZ.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-POGZ Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.