

Anti-USP25 Antibody Picoband® HRP Conjugated

Catalog Number: A06182-1-HRP

About USP25

Ubiquitin specific peptidase 25 is a protein that in humans is encoded by the USP25 gene. It is mapped to 21q21.1. Ubiquitin is a highly conserved 76-amino acid protein involved in regulation of intracellular protein breakdown, cell cycle regulation, and stress response. Ubiquitin is released from degraded proteins by disassembly of the polyubiquitin chains, which is mediated by ubiquitin-specific proteases (USPs), such as USP25.

Overview

Product Name	Anti-USP25 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Monkey, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UHP3

Technical Details

Immunogen	E.coli-derived human USP25 recombinant protein (Position: M1-D1053).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-USP25 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.