

Anti-EVI5 Antibody Picoband® Biotin Conjugated

Catalog Number: A06185-1-Biotin

About EVI5

Enables GTPase activator activity and small GTPase binding activity. Involved in positive regulation of GTPase activity and retrograde transport, endosome to Golgi. Located in cytosol.

Overview

Product Name	Anti-EVI5 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O60447

Technical Details

Immunogen	E.coli-derived human EVI5 recombinant protein (Position: N12-V810).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EVI5 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.