

Anti-EVI5 Antibody Picoband® Cy3 Conjugated

Catalog Number: A06185-1-Cy3

About EVI5

Enables GTPase activator activity and small GTPase binding activity. Involved in positive regulation of GTPase activity and retrograde transport, endosome to Golgi. Located in cytosol.

Overview

| | |
|----------------------|--|
| Product Name | Anti-EVI5 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | O60447 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human EVI5 recombinant protein (Position: N12-V810). |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-EVI5 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.