

Anti-RCC2 Antibody Picoband® PE Conjugated

Catalog Number: A06203-2-PE

About RCC2

Protein RCC2 also known as telophase disk protein of 60 kDa (TD-60) or RCC1-like protein TD-60 is a protein that in humans is encoded by the RCC2 gene. The protein encoded by this gene is a guanine exchange factor that is active on Ra1A, a small GTPase. The encoded protein and Ra1A are both essential for proper kinetochore-microtubule function in early mitosis. This protein has been shown to be a biomarker for colorectal cancer.

Overview

Product Name	Anti-RCC2 Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9P258

Technical Details

Immunogen	E.coli-derived human RCC2 recombinant protein (Position: K293-L522).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RCC2 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.