

Anti-Opioid growth factor receptor OGFR Antibody

Catalog Number: A06220


Overview

Product Name	Anti-Opioid growth factor receptor OGFR Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Opioid growth factor receptor OGFR Antibody catalog # A06220. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NZT2


Technical Details

Immunogen	Recombinant full length Human OGFR.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-Opioid growth factor receptor OGFR Antibody (A06220) Images



Immunofluorescence analysis of A549 cell using OGFR antibody. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using OGFR antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Opioid growth factor receptor OGFR Antibody

For Research Use Only. Not for use in diagnostic procedures.