

## Anti-Ogfr Antibody Picoband® APC Conjugated

Catalog Number: A06220-1-APC

### About Ogfr

Opioid growth factor receptor, also known as OGF<sub>r</sub> or the zeta-opioid receptor, is a protein which in humans is encoded by the OGFR gene. Predicted to enable G protein-coupled opioid receptor activity. Predicted to act upstream of or within regulation of growth. Predicted to be located in cytoplasm; membrane; and nucleus. Is expressed in genitourinary system and metanephros. Orthologous to human OGFR (opioid growth factor receptor).

### Overview

Product Name	Anti-Ogfr Antibody Picoband® APC Conjugated
Reactive Species	Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q99PG2

### Technical Details

Immunogen	E.coli-derived mouse Ogfr recombinant protein (Position: H217-K427).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ogfr Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.