

Anti-Ogfr Antibody Picoband® HRP Conjugated

Catalog Number: A06220-1-HRP

About Ogfr

Opioid growth factor receptor, also known as OGF_r or the zeta-opioid receptor, is a protein which in humans is encoded by the OGFR gene. Predicted to enable G protein-coupled opioid receptor activity. Predicted to act upstream of or within regulation of growth. Predicted to be located in cytoplasm; membrane; and nucleus. Is expressed in genitourinary system and metanephros. Orthologous to human OGFR (opioid growth factor receptor).

Overview

Product Name	Anti-Ogfr Antibody Picoband® HRP Conjugated
Reactive Species	Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q99PG2

Technical Details

Immunogen	E.coli-derived mouse Ogfr recombinant protein (Position: H217-K427).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ogfr Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.