

## Anti-OGFR Antibody

Catalog Number: A06220-2

### About OGFR

Opioid growth factor receptor, also known as OGFR or the zeta-opioid receptor, is a protein which in humans is encoded by the OGFR gene. Predicted to enable G protein-coupled opioid receptor activity. Predicted to act upstream of or within regulation of growth. Predicted to be located in cytoplasm; membrane; and nucleus. Is expressed in genitourinary system and metanephros. Orthologous to human OGFR (opioid growth factor receptor).

### Overview

Product Name	Anti-OGFR Antibody
Reactive Species	Human
Description	Boster Bio Anti-OGFR Antibody catalog # A06220-2. Tested in WB, IHC, ICC, IF, IP, Flow Cytometry, ELISA applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NZT2

### Technical Details

Immunogen	E.coli-derived human OGFR recombinant protein (Position: A355-P657).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 Immunoprecipitation, 1:50

Flow Cytometry (Fixed), 1-3ug/1x10<sup>6</sup> cells  
ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-OGFR Antibody

For Research Use Only. Not for use in diagnostic procedures.