

Anti-H1-2 Antibody

Catalog Number: A06243-1

Overview

| | |
|----------------------|--|
| Product Name | Anti-H1-2 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-H1-2 Antibody catalog # A06243-1. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat. |
| Application | IF, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | 500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P16403 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human H1-2 recombinant protein (Position: M1-P210). |
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-H1-2 Antibody

For Research Use Only. Not for use in diagnostic procedures.