

Anti-Leukocyte surface antigen CD53 CD53 Antibody

Catalog Number: A06249


Overview

Product Name	Anti-Leukocyte surface antigen CD53 CD53 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Leukocyte surface antigen CD53 CD53 Antibody catalog # A06249. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P19397


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CD53. AA range:91-140
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Leukocyte surface antigen CD53 CD53 Antibody (A06249) Images



Western Blot (WB) analysis of KB cells using CD53 Polyclonal antibody.



Western Blot (WB) analysis of KB cells using CD53 Polyclonal antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Leukocyte surface antigen CD53 CD53 Antibody

For Research Use Only. Not for use in diagnostic procedures.