

Anti-Fn14 TNFRSF12A Antibody

Catalog Number: A06261

About TNFRSF12A

TWEAK is a TNF family member produced by IFN-gamma-stimulated monocytes and induces multiple cell death pathways including caspase-dependent apoptosis and cathepsin B-dependent necrosis. These pathways are mediated by the binding of TWEAK by fibroblast growth factor 14 (Fn14), a member of the TNF receptor superfamily. Yeast two-hybrid experiments have shown that FN14 will bind to the TNFR-associated factors TRAF1, TRAF2, TRAF3, and TRAF5, suggesting that the signaling pathway induced by TWEAK may be mediated by one or more of these adaptor molecules. Fn14 is induced during liver regeneration and is highly expressed in hepatocellular carcinomas. In addition, TWEAK induces liver progenitor cell proliferation, suggesting Fn14 may play a role in hepatocyte growth control and liver neoplasia.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Fn14 TNFRSF12A Antibody
Reactive Species	Human
Description	Boster Bio Anti-Fn14 TNFRSF12A Antibody (Catalog # A06261). Tested in ELISA, WB, IHC-P, IF applications. This antibody reacts with Human.
Application	ELISA, IF, IHC-P, WB
Clonality	Polyclonal
Formulation	Fn14 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	Fn14 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. Avoid repeated freeze-thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Host	Rabbit
Uniprot ID	Q9NP84

Technical Details


Immunogen	Fn14 antibody was raised against a human IgG1 Fc-tagged Fn14 protein.
Predicted Reactive Species	Bovine

Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	Fn14 Antibody is affinity chromatography purified via peptide column.
Suggested Dilutions	Predicted: 62kD Observed: 68kD Optimal dilutions for each application should be determined by the researcher.


Anti-Fn14 TNFRSF12A Antibody (A06261) Images



Western blot analysis of Fn14 in HepG2 cells with Fn14 antibody at (A) 2 and (B) 4 ug/mL.



Immunohistochemistry of Fn14 in human liver tissue with Fn14 antibody at 2.5 ug/mL.



Immunofluorescence of Fn14 in Human Liver cells with Fn14 antibody at 20 ug/mL.

[Anti-Fn14 TNFRSF12A Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.