

Anti-TWEAKR TNFRSF12A Antibody

Catalog Number: A06261-2


Overview

Product Name	Anti-TWEAKR TNFRSF12A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TWEAKR TNFRSF12A Antibody catalog # A06261-2. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NP84


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

Anti-TWEAKR TNFRSF12A Antibody (A06261-2) Images



Western blot analysis of TWEAKR on different lysates using anti-TWEAKR antibody at 1/1,000 dilution. Positive control: Lane 1: NIH/3T3 Lane 2: C2C12



Immunohistochemical analysis of paraffin-embedded mouse uterus tissue using anti-TWEAKR antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TWEAKR TNFRSF12A Antibody

For Research Use Only. Not for use in diagnostic procedures.