

## Anti-BRCA1-A complex subunit RAP80 UIMC1 Antibody

Catalog Number: A06262-1


### Overview

Product Name	Anti-BRCA1-A complex subunit RAP80 UIMC1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-BRCA1-A complex subunit RAP80 UIMC1 Antibody catalog # A06262-1. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96RL1


### Technical Details

Immunogen	Recombinant fusion protein of human UIMC1(NP_057374.3).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:100

## Anti-BRCA1-A complex subunit RAP80 UIMC1 Antibody (A06262-1) Images



Western blot analysis of extracts of various cell lines, using UIMC1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunofluorescence analysis of GFP-RNF168 transgenic U2OS cells using UIMC1 antibody. Green: GFP-RNF168 fusion protein expression for DNA damage marker. Blue: DAPI for nuclear staining. RNF168 can be used to mark cells damaged by UV-A laser for they always gather around DNA damage region.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BRCA1-A complex subunit RAP80 UIMC1 Antibody

For Research Use Only. Not for use in diagnostic procedures.