

Anti-DUSP16 Antibody

Catalog Number: A06264


Overview

Product Name	Anti-DUSP16 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DUSP16 Antibody catalog # A06264. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BY84

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 577-621 of Human DUSP16
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Suggested Dilutions	WB: 1:500-1:1000

Anti-DUSP16 Antibody (A06264) Images



Western blot (WB) analysis of DUSP16 polyclonal antibody at 1:500 dilution Lane1:The Kidney tissue lysate of Rat(40ug) Lane2:The Brain tissue lysate of Mouse(40ug) Lane3:SK-OVCAR3 whole cell lysate(40ug) Lane4:A549 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DUSP16 Antibody

For Research Use Only. Not for use in diagnostic procedures.