

Anti-GPR35 Antibody

Catalog Number: A06269

About GPR35

SmarcAL1 (SWI/SNF-related matrix-associated actin-dependent regulator of chromatin a-like 1) is also known as HARP, HepA Related Protein, HHARP, SMARCA like Protein 1 and SWI/SNF Related. SmarcAL1 is a member of the SWI/SNF family of proteins. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein shows sequence similarity to the E. coli RNA polymerase-binding protein HepA. Mutations in this gene are a cause of Schimke immunosseous dysplasia (SIOD), an autosomal recessive disorder with the diagnostic features of spondyloepiphyseal dysplasia, renal dysfunction, and T-cell immunodeficiency.

Overview


Product Name	Anti-GPR35 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-GPR35 Antibody (Catalog# A06269). Tested in WB application(s). This antibody reacts with Mouse, Human.
Application	WB
Clonality	Polyclonal
Formulation	100 ug in 200 ul PBS containing 0.05% BSA and 0.05% sodium azide.
Storage Instructions	Store the antibody at 4°C, stable for 6 months. For long-term storage, store at -20°C. Avoid repeated freeze and thaw cycles.
Host	Rabbit
Uniprot ID	Q9HC97

Technical Details


Immunogen	A partial length recombinant GPR35 protein (amino acids 173-305)
Predicted Reactive Species	Chimpanzee
Cross Reactivity	No cross reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Purification	Protein A Chromatography
Suggested Dilutions	Western blot: 2-4 µg/ml

Anti-GPR35 Antibody (A06269) Images



Western blot analysis of GPR35. Anti-GPR35 antibody was used at 2 ug/ml on HCT-116 lysate.



Western blot analysis of GPR35. Anti-GPR35 antibody was used at 1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GPR35 Antibody