

## Anti-Zinc finger protein 521 ZNF521 Antibody

Catalog Number: A06292

### About ZNF521

The zinc finger protein 521 (ZNF521) is a transcription factor containing an N-terminal transcriptional repressor motif and 30 zinc finger domains. It plays a role in both erythroid cell and osteoblast differentiation during development, inhibiting the activities of early B-cell factor 1 (EBF1) in erythroid cells and Runx2 in osteoblast precursors. ZNF521 binds to both Runx2 and histone deacetylase 3 (HDAC3), promotes their association and antagonizes Runx2 transcriptional activity in a HDAC3-dependent manner, thereby regulating osteoblast differentiation, skeletal development, and bone homeostasis.

### Overview

Product Name	Anti-Zinc finger protein 521 ZNF521 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Zinc finger protein 521 ZNF521 Antibody (Catalog # A06292). Tested in ELISA, IHC-P applications. This antibody reacts with Human.
Application	ELISA, IHC-P
Clonality	Polyclonal
Formulation	ZNF521 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	ZNF521 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. Avoid repeated freeze-thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Host	Rabbit
Uniprot ID	Q96K83

### Technical Details


Immunogen	ZNF521 antibody was raised against a 19 amino acid synthetic peptide near the carboxy terminus of human ZNF521. The immunogen is located within amino acids 1160 - 1210 of ZNF521.
Predicted Reactive Species	Mouse
Cross Reactivity	ZNF521 antibody is human specific.
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	ZNF521 Antibody is affinity chromatography purified via peptide column.

Suggested Dilutions

ZNF521 antibody can be used for detection of ZNF521 by immunohistochemistry at 5 ug/mL.  
Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.

## Anti-Zinc finger protein 521 ZNF521 Antibody (A06292) Images

---



Immunohistochemistry of ZNF521 in human kidney tissue with ZNF521 antibody at 5 ug/mL.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zinc finger protein 521 ZNF521 Antibody

For Research Use Only. Not for use in diagnostic procedures.