

## Anti-Protein NOV homolog NOV Antibody

Catalog Number: A06319

### Overview


Product Name	Anti-Protein NOV homolog NOV Antibody
Reactive Species	Human
Description	Boster Bio Anti-Protein NOV homolog NOV Antibody catalog # A06319. Tested in IHC applications. This antibody reacts with Human.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P48745

### Technical Details

Immunogen	Recombinant fusion protein of human NOV(NP_002505.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:100

## Anti-Protein NOV homolog NOV Antibody (A06319) Images

---



Immunohistochemistry of paraffin-embedded human prostate using NOV antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein NOV homolog NOV Antibody

For Research Use Only. Not for use in diagnostic procedures.