

Anti-K-Cadherin-6 CDH6-Antibody

Catalog Number: A06353


Overview

Product Name	Anti-K-Cadherin-6 CDH6-Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-K-Cadherin-6 CDH6-Antibody catalog # A06353. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P55285

Technical Details

Immunogen	Synthesized peptide derived from K-cadherin . at AA range: 180-260
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using a epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-300 ELISA 1:20000 IF 1:50-200

Anti-K-Cadherin-6 CDH6-Antibody (A06353) Images



Western Blot analysis of BT474, A549, SHSY5Y, 293 cells using K-cadherin Polyclonal Antibody. Secondary antibody was diluted at 1:20000

1 Publications Citing This Product

1. PubMed ID: 10.1002/1878-0261.12947, CDH6 β activated α 1 β 3 crosstalks with α 2 β 1 to trigger cellular adhesion and invasion in metastatic ovarian and renal cancers

Visit bosterbio.com/anti-k-cadherin-antibody-a06353-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-K-Cadherin-6 CDH6-Antibody

For Research Use Only. Not for use in diagnostic procedures.