

## Anti-DGK-Alpha DGKA Antibody

Catalog Number: A06354-2


### Overview

Product Name	Anti-DGK-Alpha DGKA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DGK-Alpha DGKA Antibody catalog # A06354-2. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23743

### Technical Details

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DGKA (NP_958852.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:100-1:500

## Anti-DGK-Alpha DGKA Antibody (A06354-2) Images



Western blot analysis of extracts of various cell lines, using DGKA pAb antibody at 1:500 dilution.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DGK-Alpha DGKA Antibody

For Research Use Only. Not for use in diagnostic procedures.