

Anti-MT2-MMP/MMP15 Antibody Picoband® HRP Conjugated

Catalog Number: A06396-1-HRP

About MMP15

Matrix metalloproteinase 15 also known as MMP15 is an enzyme that in humans is encoded by the MMP15 gene. This gene encodes a member of the peptidase M10 family and membrane-type subfamily of matrix metalloproteinases (MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Members of this subfamily contain a transmembrane domain suggesting that these proteins are expressed at the cell surface rather than secreted. The encoded preproprotein is proteolytically processed to generate the mature protease. This protein may play a role in cancer progression.

Overview

| | |
|----------------------|---|
| Product Name | Anti-MT2-MMP/MMP15 Antibody Picoband® HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P51511 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human MT2-MMP/MMP15 recombinant protein (Position: Y132-A202). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MT2-MMP/MMP15 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.