

## Anti-VPS28 Antibody

Catalog Number: A06407-2

### About VPS28

VPS28 is a subunit of the ESCRT-I complex (endosomal sorting complex I required for transport), which functions in the transport and sorting of proteins into subcellular vesicles. The ESCRT-I complex is composed of VPS23 (also known as TSG101), VPS28, VPS37, and MVB12A or MVB12B. VPS28 interacts with VPS23. Two alternative transcripts encoding different VPS28 isoforms have been described.

### Overview

Product Name	Anti-VPS28 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-VPS28 Antibody catalog # A06407-2. Tested in WB, ICC/IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UK41

### Technical Details

Immunogen	E.coli-derived human VPS28 recombinant protein (Position: H3-F218).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-VPS28 Antibody

For Research Use Only. Not for use in diagnostic procedures.