

Anti-SIN3B Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A06424-2-Fluoro488

About SIN3B

Predicted to enable transcription corepressor activity. Predicted to be involved in cardiac muscle tissue development; negative regulation of transcription by RNA polymerase II; and skeletal muscle tissue development. Predicted to be located in nucleus. Predicted to be part of Sin3-type complex. Predicted to be active in chromatin.

Overview

Product Name	Anti-SIN3B Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O75182

Technical Details

Immunogen	E.coli-derived human SIN3B recombinant protein (Position: Q127-E1127). Human SIN3B shares 86.5% amino acid (aa) sequence identity with mouse SIN3B.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SIN3B Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.