

Anti-SIN3B Antibody Picoband® HRP Conjugated

Catalog Number: A06424-2-HRP

About SIN3B

Predicted to enable transcription corepressor activity. Predicted to be involved in cardiac muscle tissue development; negative regulation of transcription by RNA polymerase II; and skeletal muscle tissue development. Predicted to be located in nucleus. Predicted to be part of Sin3-type complex. Predicted to be active in chromatin.

Overview

Product Name	Anti-SIN3B Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O75182

Technical Details

Immunogen	E.coli-derived human SIN3B recombinant protein (Position: Q127-E1127). Human SIN3B shares 86.5% amino acid (aa) sequence identity with mouse SIN3B.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SIN3B Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.