

Anti-Fc receptor-like protein 5 FCRL5 Antibody

Catalog Number: A06435-1


Overview

Product Name	Anti-Fc receptor-like protein 5 FCRL5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Fc receptor-like protein 5 FCRL5 Antibody catalog # A06435-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96RD9


Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human FCRL5.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Fc receptor-like protein 5 FCRL5 Antibody (A06435-1) Images



Western blot analysis of CD307e expression in Mouse Liver (A) whole cell lysates.



Western blot analysis of CD307e expression in Jurkat (A), NS-1 (B), PC12 (C) whole cell lysates.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fc receptor-like protein 5 FCRL5 Antibody

For Research Use Only. Not for use in diagnostic procedures.