

Anti-CD218b IL18RAP Antibody

Catalog Number: A06457


Overview

Product Name	Anti-CD218b IL18RAP Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CD218b IL18RAP Antibody catalog # A06457. Tested in IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95256


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 291-340
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

Anti-CD218b IL18RAP Antibody (A06457) Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

1 Publications Citing This Product

1. PubMed ID: 10.1016/j.exer.2021.108750, Wedelolactone ameliorates Pseudomonas aeruginosa-induced inflammation and corneal injury by suppressing caspase-4/5/11/GSDMD-mediated non-canonical pyroptosis

Visit bosterbio.com/anti-cd218b-antibody-a06457-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD218b IL18RAP Antibody

For Research Use Only. Not for use in diagnostic procedures.