

## Anti-Legumain LGMN Antibody

Catalog Number: A06476


### Overview

Product Name	Anti-Legumain LGMN Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Legumain LGMN Antibody catalog # A06476. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99538

### Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 60-140
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	LGMN antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

## Anti-Legumain LGMN Antibody (A06476) Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

## 1 Publications Citing This Product

1. PubMed ID: 32229193, Chen HC,Rui W,You SY,Liu XW,Huang J,Chen HY.Evaluation of the anti-cervical cancer effect of a prodrug 7CBZ-AAN-DOX with hypoxic cell culture and tumor-bearing zebrafish models. Exp Cell Res.2020 Jun 1;391(1):111980.doi: 10.1016/j.yexcr.2020.111980.Epub 2

Visit [bosterbio.com/anti-legumain-lgm-antibody-a06476-boster.html](https://bosterbio.com/anti-legumain-lgm-antibody-a06476-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Legumain LGMN Antibody

For Research Use Only. Not for use in diagnostic procedures.