

Anti-CSRP3 Antibody

Catalog Number: A06522-2

About CSRP3

CSRP3/MLP (cysteine-rich protein 3/muscle Lim protein), a member of the cysteine-rich protein family, is a muscle-specific LIM-only factor specifically expressed in skeletal muscle. The cytoplasmic CSRP3 serves as a scaffold protein that primarily interacts and co-localizes with different Z-disc components of myofibrils (PMID: 31979369, 17519359).

Overview

Product Name	Anti-CSRP3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CSRP3 Antibody catalog # A06522-2. Tested in WB, IHC, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P50461

Technical Details

Immunogen	E.coli-derived human CSRP3 recombinant protein (Position: N3-E191).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CSR3 Antibody

For Research Use Only. Not for use in diagnostic procedures.