

## Anti-Transaldolase TALDO1 Antibody

Catalog Number: A06530


### Overview

Product Name	Anti-Transaldolase TALDO1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Transaldolase TALDO1 Antibody catalog # A06530. Tested in WB,IF applications. This antibody reacts with Human,Mouse.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P37837

### Technical Details

Immunogen	Recombinant full length Human TALDO1.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

## Anti-Transaldolase TALDO1 Antibody (A06530) Images



WesternBlot (WB) analysis of TALDO1 polyclonal antibody

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transaldolase TALDO1 Antibody

For Research Use Only. Not for use in diagnostic procedures.