

Anti-Urotensin II UTS2 Antibody

Catalog Number: A06633

Overview

Product Name	Anti-Urotensin II UTS2 Antibody
Reactive Species	Human, Monkey
Description	Boster Bio Anti-Urotensin II UTS2 Antibody catalog # A06633. Tested in WB, ELISA applications. This antibody reacts with Human, Monkey.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95399

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human Urotensin II. AA range:48-97
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-Urotensin II UTS2 Antibody (A06633) Images



Western Blot analysis of various cells using Urotensin II Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Western blot analysis of lysate from 293 cells, using Urotensin II antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Urotensin II UTS2 Antibody

For Research Use Only. Not for use in diagnostic procedures.